

# TAKE CONTROL OF YOUR IHC STAINING



Cell Marque™ offers an extensive array of positive control slides that correspond with the majority of our antibody portfolio. Each tissue block is selected and reviewed by a certified pathologist in order to provide our customers with quality assurance and the best overall experience when testing and validating their antibodies.

## The Cell Marque™ control slide portfolio offers:

- FFPE human tissues
- Positively charged slides
- Tissue cut optimally at 4μ
- Guaranteed antigenicity
- 36 month shelf life
- 5 unstained slides per pack
- Availability of rare and viral infected tissues
- Increased confidence in antibody performance
- Assistance with validation of new antibodies
- Improved quality control of internal tissue



©2015 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH and Cell Marque are trademarks of Sigma-Aldrich Co. LLC.