





Clone Wars: SOX-10

mouse monoclonal vs. rabbit monoclonal, clone: EP268

Top-left: mouse monoclonal SOX-10 on desmoplastic melanoma 400X

Top-right: rabbit monoclonal SOX-10 (EP268) on desmoplastic melanoma 400X



Bottom-left: mouse monoclonal SOX-10 on desmoplastic melanoma 200X

Bottom-right: rabbit monoclonal SOX-10 (EP268) on desmoplastic melanoma 200X

SOX-10 (EP268)

Volume.....	Part No.
0.1 ml, concentrate.....	383R-14
0.5 ml, concentrate.....	383R-15
1 ml, concentrate.....	383R-16
1 ml, prediluted.....	383R-17
7 ml, prediluted.....	383R-18
Positive control slides ..	3835

Benefits of Rabbit Monoclonal SOX-10:

- *In vitro* diagnostic
- Nuclear visualization to eliminate cytoplasmic melanin interference
- Sensitivity of a rabbit antibody with the specificity and cleanliness of a monoclonal antibody for a strong signal to noise ratio
- More intense staining than mouse counterparts
- Labels desmoplastic melanoma with a higher sensitivity and specificity than anti-S-100
- Rabbit antibody is compatible with standard automation and detection used in diagnostic IHC

CELL MARQUE

RabMAb®
Technology from Abcam

Cell Marque™ – a Sigma-Aldrich® Company

Toll-Free: 800.665.7284
Direct: 916.746.8900

service@cellmarque.com
www.cellmarque.com