

30x Faster!

Drastically reduce staining time compared to standard procedures with the new Fontana-Masson Stain Kit.

- Microwave procedure reduces melanin incubation time from 30 minutes to 1 minute.
- Argentaffin staining incubation time is decreased from 60 minutes to 2 minutes.
- Ammoniacal Silver Solution is ready to use immediately upon preparation - no lengthy wait time.
- Each component is performance tested to provide quality products you can trust.



Below left: Melanin stain.

Below right: Argentaffin stain

Fontana-Masson Stain Kit Cat. No. HT200

Fontana Masson Stain Kit Components:
Silver Nitrate Solution (Cat. No. HT2001- 50 ml),
Gold Chloride Solution (Cat. No. HT2002 - 125 ml),
Sodium Thiosulfate Solution (Cat. No. HT2003 - 125 ml)

service@cellmarque.com • www.sigmaaldrich.com



The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

©2017 Sigma-Aldrich Co. LLC. All rights reserved. Sigma-Aldrich Co., LLC is a subsidiary of Merck KGaA, Darmstadt, Germany. MilliporeSigma and the Vibrant M are trademarks of Merck KGaA, Darmstadt, Germany. Sigma-Aldrich is a trademark of Sigma-Aldrich Co. LLC or its affiliates.
Proj. No. 17-05-0114-PRT-FLYER Ver: 0.0 05/2017

**MILLIPORE
SIGMA**